

# ANPz102/UHV - linear z-nanopositioner

1009457

## Technical Specifications

<b>Size and Dimensions</b>	
footprint; height	24 mm x 24 mm; 27 mm
maximum installation space	24 mm x 24 mm; 32 mm
weight	53 g
height	27 mm
<b>Materials</b>	
positioner body	titanium
actuator	PZT ceramics
connecting wires	insulated twisted pair, copper
<b>Coarse Positioning Mode</b>	
travel range (step mode)	4.8 mm
maximum drive velocity @ 300 K	~ 3 mm/s
input voltage range	0 - 60 V
<b>Fine Positioning Mode</b>	
fine linear positioning range @ 300 K	3.5 $\mu$ m
fine positioning resolution	sub-nm
input DC voltage range @ 300 K	0 - 100 V
<b>Load (@ ambient conditions)</b>	
maximum load	200 g
maximum dynamic force along the axis	5 N
<b>General Specifications</b>	
environment	/UHV

